



officine
orobiche

quantity

1

off No.

P19/4121 pos 1 of 1 rev 0

model

CV25.MO.MO.M14.WP.XX

tag

s1
s2
s3
s4
s5

s6
s7
s8
s9
s10

s11
s12
s13
s14
s15

s16
s17
s18
s19
s20

01
02
03
04
05

body mat. MONEL
connections 1" ANSI 300 RF 125-250 AARH
float MONEL
attraction sleev Nickel
flow indicator NO

21
22
23
24
25

body rating
max-press.
max-temp.
fluid

10 barg
100 °C
SEA WATER

06
07
08

housing ALUMINIUM WP IP 66
AREA .

09
10
11
12

electric connection ISO M20 x 1,5
switch groups 1
microswitch M14 (inert gas hermetically sealed)
contact SPDT (48 V DC 2 A)

26
27
28
29
30
31
32

nominal flow
range
set (min)
diff.
normal flow
adjust.(min)
scale

3 m3/h
1---2.5 m3/h
0,9 m3/h
0,9 m3/h
2,4 m3/h
supplier
.

13
14
15
16
17

connector
cable gland
fins
label
paint

AISI 316

33
34
35
36
37

nace

NO

18
19
20

tube MONEL
internal parts MONEL
antisteam

38
39
40

mat_cone


MONEL

v1
v2
v3
v4
v5


n1
n2
n3
n4

NORMAL FLOW RATE: 40 L/MIN.
SET FOR LOW FLOW CALIBRATED AT 15 L/MIN
F.70580
modificato tag con F.74581

X



.



mag-19

issued UV.....date

approved UT.....date